

Mandatory setup for WBC Differential Counter for Windows XP and Windows 7

If you use WBC Differential Counter on Windows XP or Windows 7 you must do a mandatory onetime setup. Before doing this setup the program probably looks unreadable as shown on Figure 1.




Figure 1.

After you downloaded and ran the program, click over the icon with the right mouse button to open Properties – Figure 2.




Figure 2.

After that on the new window choose Layout . Enter on the Screen Buffer Size Width : 120 and enter on Screen size Width : 120 and click OK – Figure 3. After that you have to restart the program.




Figure 3.

After you have done the setting the program should look fine as in Figure 4.




Figure 4.